


## Aero Spacer (Valved Holding Chamber)

Model: AS175



- High efficiency one way valve for easy inhalation and exhalation
- Universal Backpiece, can be used with all commonly prescribed Metered Dose Inhalers
- Anti-static spacer improves the delivery of MDI by inhibiting the adherence of drug to spacer wall
- Soft and flexible silicone mask provides comfortable touch and reliable seal
- Portable design & easy holding
- Easy disassembled and reassembled for cleaning
- Mouthpiece cap can protect the mouthpiece from dirt and debris when the device is not in use
- This product contains no latex and BPA free.
- Small, medium & large facemask to cater to the different needs.
- Small facemask: 0-10 months
- Medium facemask: 12 months – 5 yrs.
- Large facemask or mouthpiece: 5 yrs. +
- Chamber volume: 175ml
- For single patient use only

Art No	Qty/Box	Country	CPA/CE	SPF1	40°C1	40°C2
with Mouthpiece & Cap (3 yrs. +)	50	10.000	0.130	19.620	10.130	49.830
with Small Facemask (0-10 months)	50	11.000	0.130	19.620	10.130	49.830
with Medium Facemask (12 months - 5 yrs.)	50	11.000	0.130	19.620	10.130	49.830
with Large Facemask (5 yrs. +)	50	11.000	0.130	19.620	10.130	49.830

 **Dove and Temperature Control**  
 20°C to 40°C with an extra  
 40°C/104°F version  
[www.rossmax.com](http://www.rossmax.com)  
 0800 828 828 ext. 2000